

CareSens™ **N** **Premier** Blood Glucose Monitoring System

- 2.5" Large, Backlit Display
- Efficient Diabetes Management with Tagging Symbols
- Strip Ejector for Easy Ejection of Used Test Strips



Product specifications

Measurement range	20 – 600 mg/dL (1.1–33.3 mmol/L)
Sample size	Minimum 0.5 µL
Test time	5 seconds
Sample type	Fresh capillary whole blood
Calibration	Plasma-equivalent
Assay method	Electrochemical
Battery life	3,000 tests
Power	Two 3.0 V lithium batteries (disposable, type CR2032)
Memory	1,000 test results
Size	103 x 54 x 17 (mm)
Weight	70.9 g (with batteries)

